

Product Datasheet

Noggin (Noggin) antibody (orb1143076)

Description	Mouse monoclonal antibody to Noggin (Noggin)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 28-232 of human noggin (NP_005441) produced in E.coli.
Preservatives	PBS (pH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	23.7 kDa(Predicted)
Uniprot ID	Q13253
Dilution Range	WB 1:20000~40000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.