

Product Datasheet

MAPK13 (SAPK4) antibody (orb1143386)

Description	Mouse monoclonal antibody to MAPK13 (SAPK4)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human MAPK13(NP_002745) produced in HEK293T cell.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	41.9 kDa(Predicted)
Uniprot ID	O15264
Dilution Range	WB 1:200~2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.