

## Product Datasheet

### HSPA8 antibody (orb1143694)

<b>Description</b>	Mouse monoclonal antibody to HSPA8
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	Full length human recombinant protein of human HSPA8(NP_006588) produced in HEK293T cell.
<b>Preservatives</b>	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>MW</b>	70.7 kDa(Predicted)
<b>Uniprot ID</b>	<b>V9HW22</b>
<b>Dilution Range</b>	WB 1:1000, IHC 1:150, IF 1:100
<b>Expiration Date</b>	12 months from date of receipt.