

Product Datasheet

HOMER1 antibody (orb1143734)

Description	Mouse monoclonal antibody to HOMER1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-354 of human HOMER1 (NP_004263) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	40.1 kDa(Predicted)
Uniprot ID	Q5U5K4
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.