

Product Datasheet

GTF2B antibody (orb1143824)

Description	Mouse monoclonal antibody to GTF2B
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 43-316 of human GTF2B(NP_001505) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	34.7 kDa(Predicted)
Uniprot ID	Q00403
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.