

Product Datasheet

GPHN (Gephyrin) antibody (orb1143856)

Description	Mouse monoclonal antibody to GPHN (Gephyrin)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human GPHN (NP_065857) produced in HEK293T cell.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	83.3 kDa(Predicted)
Uniprot ID	Q9NQX3
Dilution Range	WB 1:1000~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.