

Product Datasheet

FGFR2 antibody (orb1143999)

Description	Mouse monoclonal antibody to FGFR2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human FGFR2 (NP_000132) produced in HEK293T cell.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	89.7 kDa(Predicted)
Uniprot ID	P21802
Dilution Range	WB 1:200 - 1:1000, FLOW 1:100
Expiration Date	12 months from date of receipt.