

Product Datasheet

EXOSC3 antibody (orb1144059)

Description	Mouse monoclonal antibody to EXOSC3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human EXOSC3(NP_057126) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	29.4 kDa(Predicted)
Uniprot ID	Q9NQT5
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.