

Product Datasheet

CTNND1 (Delta 1 Catenin) antibody (orb1144347)

Description	Mouse monoclonal antibody to CTNND1 (Delta 1 Catenin)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human CTNND1 (NP_001078927) produced in HEK293T cell.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	108 kDa(Predicted)
Uniprot ID	O60716
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.