

Product Datasheet

c-Myc (c-Myc) antibody (orb1144434)

Description	Mouse monoclonal antibody to c-Myc (c-Myc)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human Myc (NP_002458) produced in E.coli.
Preservatives	PBS (pH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	50.4 kDa(Predicted)
Uniprot ID	P01106
Dilution Range	WB 1:50~100, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.