

Product Datasheet

CCNB1IP1 antibody (orb1144570)

Description	Mouse monoclonal antibody to CCNB1IP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human CCNB1IP1 (NP_878269) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	31.4 kDa(Predicted)
Uniprot ID	Q9NPC3
Dilution Range	WB 1:500, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.