

Product Datasheet

ARHGDI (RhoGDI) antibody (orb1144799)

Description	Mouse monoclonal antibody to ARHGDI (RhoGDI)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human ARHGDI (NP_004300) produced in HEK293T cell.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	23.0 kDa(Predicted)
Uniprot ID	V9HWE8
Dilution Range	WB 1:200~1000, FLOW 1:100
Expiration Date	12 months from date of receipt.