

Product Datasheet

AP2M1 antibody (orb1144833)

Description	Mouse monoclonal antibody to AP2M1
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 97-383 of human AP2M1(NP_001020376) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	49.2 kDa(Predicted)
Uniprot ID	B4DNB9
Dilution Range	WB 1:200~500, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.