

Product Datasheet

AP2M1 antibody (orb1144835)

Description	Mouse monoclonal antibody to AP2M1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human AP2M1 (NP_001020376) produced in HEK293T cell.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	49.2 kDa(Predicted)
Uniprot ID	B4DNB9
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.