

Product Datasheet

SLC18A2 (VMAT2) antibody (orb1144884)

Description	Mouse monoclonal antibody to SLC18A2 (VMAT2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human SLC18A2 expression vector
Preservatives	PBS (pH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	55.7 kDa(Predicted)
Uniprot ID	Q05940
Dilution Range	WB 1:500 IHC 1:50
Expiration Date	12 months from date of receipt.