

Product Datasheet

SERPINA1 antibody (orb1144891)

Description	Mouse monoclonal antibody to SERPINA1
Species/Host	Mouse
Reactivity	Canine, Human, Monkey
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human SERPINA1 expression vector
Preservatives	PBS (pH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	46.7 kDa(Predicted)
Uniprot ID	P01009
Dilution Range	WB 1:500~1000, IHC 1:150, FLOW 1:100, IP 2ug/500ul
Expiration Date	12 months from date of receipt.