

Product Datasheet

PLK1 antibody (orb1144928)

Description	Mouse monoclonal antibody to PLK1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human PLK1 expression vector
Preservatives	PBS (pH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	68.1 kDa(Predicted)
Uniprot ID	P53350
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100, IP 2ug/500ul
Expiration Date	12 months from date of receipt.