

Product Datasheet

NT3 (Neurotrophin 3) antibody (orb1144958)

Description	Mouse monoclonal antibody to NT3 (Neurotrophin 3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 139-257 of human NT3
Preservatives	PBS (pH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	27.4 kDa(Predicted)
Uniprot ID	P20783
Dilution Range	WB: 1:5000, IHC 1:50
Expiration Date	12 months from date of receipt.