

Product Datasheet

ARHGDIA (RhoGDI) antibody (orb1145118)

Description	Mouse monoclonal antibody to ARHGDIA (RhoGDI)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human ARHGDIA(NP_004300) produced in HEK293T cell.
Preservatives	PBS (pH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	23 kDa(Predicted)
Uniprot ID	V9HWE8
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.