

Product Datasheet

PKD1 (phospho-S910) antibody [Out of stock] (orb1145396)

Description	Rabbit polyclonal antibody to PKD1 (phospho-S910)
Species/Host	Rabbit
Tested Applications	IHC, WB
Immunogen	Synthetic phosphopeptide derived from human PKD1 around the phosphorylation site of Serine 910.
Preservatives	1 mg/ml in Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.2.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
MW	~ 115 kDa
Uniprot ID	#N/A
Dilution Range	WB: 1:500~1:1000 IHC: 1:50~1:200
Expiration Date	12 months from date of receipt.