

Product Datasheet

JNK1/2/3 (Phospho-T183/Y185) antibody (orb1145692)

Description	Rabbit polyclonal antibody to JNK1/2/3 (Phospho-T183/Y185)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic phosphopeptide derived from human JNK1/2/3 around the phosphorylation site of Threonine 183 / Tyrosine 185.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 54 kDa
Uniprot ID	P45983
Dilution Range	WB: 1:10000~1:20000 IHC: 1:50~1:200 IF: 1:50~1:200
Expiration Date	12 months from date of receipt.