

Product Datasheet

ENPP1/PC1 Antibody (orb1145755)

Description

ENPP1/PC1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, FC, IHC, WB

Immunogen

E.coli-derived human ENPP1/PC1 recombinant protein (Position: E347-N789).

Preservatives

4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na₂HPO₄

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

130-150 kDa

Uniprot ID


P22413

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human Immunohistochemistry(Paraffin-embedded Section), 1-2 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human Direct ELISA, 0.1-0.5 µg/ml, Human

Expiration Date


12 months from date of receipt.



IHC analysis of ENPP1/PC1 using anti-ENPP1/PC1 antibody



IHC analysis of ENPP1/PC1 using anti-ENPP1/PC1 antibody



IHC analysis of ENPP1/PC1 using anti-ENPP1/PC1 antibody