


Product Datasheet

Poliovirus Receptor/PVR Antibody (monoclonal, 5I13D1) (orb1145772)

Description

Poliovirus Receptor/PVR Antibody (monoclonal,


Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC, WB

Applications

IHC analysis
of Poliovirus
Receptor/PVR

Immunogen

E.coli-derived human Poliovirus Receptor/PVR recombinant protein (Position: D28-E331).

...

Form/Appearance


Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.



IHC analysis
of Poliovirus
Receptor/PVR


Note

For research use only

...

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.



IHC analysis
of Poliovirus
Receptor/PVR

Isotype

Mouse IgG2a

Clonality

Monoclonal

MW

70-80 kDa

Uniprot ID

P15151

...

Dilution Range

Western blot, 0.25-0.5 µg/ml, Human Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Mouse Flow Cytometry, 1-3 µg/1x10⁶ cells, Human

Expiration Date

12 months from date of receipt.