

Product Datasheet

PC4/SUB1 Antibody (monoclonal, 6B5B10) (orb1145777)

Description

PC4/SUB1 Antibody (monoclonal, 6B5B10)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

E.coli-derived human PC4/SUB1 recombinant protein (Position: N62-L127).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW

18 kDa

Uniprot ID



[P53999](#)

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Mouse, Rat

Expiration Date


12 months from date of receipt.



IHC analysis of PC4/SUB1 using anti-PC4/...



IHC analysis of PC4/SUB1 using anti-PC4/...



IHC analysis of PC4/SUB1 using anti-PC4/...