

Product Datasheet

Cdk2 Antibody (monoclonal, 5B12D1) (orb1145854)

Description

Cdk2 Antibody (monoclonal, 5B12D1)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, IHC, WB

Applications
Immunogen

E.coli-derived human Cdk2 recombinant protein (Position: E81-L298). Human Cdk2 shares 98.6% amino acid (aa) sequence identity with rat Cdk2.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2a

Clonality

Monoclonal

MW

30 kDa

Uniprot ID


[P24941](#)

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Mouse

Expiration Date


12 months from date of receipt.



IHC analysis of Cdk2 using anti-Cdk2 ant...



IHC analysis of Cdk2 using anti-Cdk2 ant...



IHC analysis of Cdk2 using anti-Cdk2 ant...