

Product Datasheet

MSH2 Antibody (monoclonal, 6B4F7) (orb1145866)

Description

MSH2 Antibody (monoclonal, 6B4F7)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ICC, IF, IHC, WB

Applications

Immunogen

E.coli-derived human MSH2 recombinant protein (Position: Q337-N583). Human MSH2 shares 94% and 93% amino acid (aa) sequence identity with mouse and rat MSH2, respectively.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW

105 kDa

Uniprot ID


[P43246](#)

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human
Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human
Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human

Expiration Date


12 months from date of receipt.



IHC analysis of MSH2 using anti-MSH2 ant...



IHC analysis of MSH2 using anti-MSH2 ant...



IHC analysis of MSH2 using anti-MSH2 ant...