

Product Datasheet

ACE2 antibody (Internal) (orb1146446)

Description	ACE2 antibody (Internal)
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ACE2 (Internal)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 100kDa band was also observed on full-length recombinant Human ACE2 protein (calculated MW of 92.4kDa according to NP_001358344.1). Recommended concentration 1ug/ml. Primary incubation 1 hour at room temperature.
Protein Sequence	PROPRIETARY;
Clonality	Polyclonal
MW	92.4
Entrez	59272
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000.
Expiration Date	12 months from date of receipt.