

Product Datasheet

Mouse TRACP-5b(Tartrate Resistant Acid Phosphatase 5b) ELISA Kit (orb1146716)

Description The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

Reactivity Mouse

Range 0.78-50 ng/mL

Note For research use only

Sample Types serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Sensitivity 0.47 ng/mL

Expiration Date Please enquire.

