

Product Datasheet

Mouse VAP-1(VascularAdhesion Protein 1) ELISA Kit (orb1146724)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	O70423
Sensitivity	0.10 ng/mL
Expiration Date	Please enquire.

