

Product Datasheet

Mouse CXCL15(Chemokine C-X-C-Motif Ligand 15) ELISA Kit (orb1146727)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	78.13-5000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9WVL7
Sensitivity	46.88 pg/mL
Expiration Date	Please enquire.

