

Product Datasheet

Human STUB1(E3 ubiquitin-protein ligase CHIP) ELISA Kit (orb1146733)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	31.25-2000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9UNE7
Sensitivity	18.75 pg/mL
Expiration Date	Please enquire.

