

Product Datasheet

Human mTOR(MammalianTarget of Rapamycin) ELISA Kit (orb1146741)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P42345
Sensitivity	0.10 ng/mL
Expiration Date	Please enquire.

