

Product Datasheet

Human P-Cadherin(Cadherin, Placental) ELISA Kit (orb1146744)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	62.50-4000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P22223
Sensitivity	37.50 pg/mL
Expiration Date	Please enquire.

