

## Product Datasheet

# Mouse MDGA1(MAM domain-containing glycosylphosphatidylinositol anchor protein 1) ELISA Kit (orb1146774)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>      | Mouse                                                                                                   |
| <b>Range</b>           | 0.16-10 ng/mL                                                                                           |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Sample Types</b>    | tissue homogenates, cell lysates and other biological fluids                                            |
| <b>Uniprot ID</b>      | <b>Q0PMG2</b>                                                                                           |
| <b>Sensitivity</b>     | 0.054 ng/mL                                                                                             |
| <b>Expiration Date</b> | Please enquire.                                                                                         |

