

Product Datasheet

Mouse BCAT2(Branched-chain-amino-acid aminotransferase, mitochondrial) ELISA Kit (orb1146805)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	15.63-1000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	035855
Sensitivity	6.5 pg/mL
Expiration Date	Please enquire.

