

Product Datasheet

Mouse SIRT2(Sirtuin 2) ELISA Kit (orb1146825)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	1.57-100 ng/mL
Note	For research use only
Sample Types	tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q8VDQ8
Sensitivity	0.67 ng/mL
Expiration Date	Please enquire.

