

Product Datasheet

Human NEDD4-L(E3 ubiquitin-protein ligase NEDD4-like) ELISA Kit (orb1146830)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.32-20 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids
Uniprot ID	Q96PU5
Sensitivity	0.119 ng/mL
Expiration Date	Please enquire.

