

Product Datasheet

Rat TRPV1(Transient Receptor Potential Cation Channel Subfamily V, Member 1) ELISA Kit (orb1146838)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	Tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	035433
Sensitivity	0.059 ng/mL
Expiration Date	Please enquire.

