

Product Datasheet

Mouse STAT3(Signal Transducer And Activator Of Transcription 3) ELISA Kit (orb1146854)

Description The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

Reactivity Mouse

Range 0.16-10 ng/mL

Note For research use only

Sample Types tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Uniprot ID [P42227](#)

Sensitivity 0.062 ng/mL

Expiration Date Please enquire.

