

Product Datasheet

Mouse CDNF(Cerebral Dopamine Neurotrophic Factor) ELISA Kit (orb1146917)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	31.25-2000 pg/mL
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8CC36
Sensitivity	15.3 pg/mL
Expiration Date	Please enquire.

