

Product Datasheet

Rat STAT3(Signal Transducer And Activator Of Transcription 3) ELISA Kit (orb1146919)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P52631
Sensitivity	0.062 ng/mL
Expiration Date	Please enquire.

