

Product Datasheet

Hamster IFN γ (Interferon Gamma) ELISA Kit (orb1146926)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Hamster
Range	15.63-1000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	O35497
Sensitivity	5.9 pg/mL
Expiration Date	Please enquire.

