

Product Datasheet

Rat CDKN2A(Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit (orb1146952)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.32-20 ng/mL
Note	For research use only
Sample Types	Tissue homogenates and other biological fluids.
Uniprot ID	Q9R0Z3
Sensitivity	0.108 ng/mL
Expiration Date	Please enquire.

