

Product Datasheet

Human OLIG2(Oligodendrocyte Lineage Transcription Factor 2) ELISA Kit (orb1146965)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	15.63-1000 pg/mL
Note	For research use only
Sample Types	Tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q13516
Sensitivity	5.7 pg/mL
Expiration Date	Please enquire.

