

Product Datasheet


Human IGFBP-3 Protein (orb1146966)

Description

Human IGFBP-3, His Tag (IG3-H52H9) is expressed from

Human IGFBP-3, His Tag ELISA
0.2 µg of Biotinylated human IGF-II, Anti-IGF-3 Ig per well

Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IGFBP-3
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 30.6 kDa. The protein migrates as 40-50 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Source	Human IGFBP-3, His Tag (orb1146966) is expressed from human 293 cells (HEK293). It contains AA Gly 28 - Lys 291 (Accession # P17936-1).
MW	30.6 kDa
Uniprot ID	P17936
NCBI	NP_000589.2
Expiration Date	6 months from date of receipt.



Immobilized
Biotinylated
Human IGF-
II, A...