

Product Datasheet

Human PDGF R beta / CD140b Protein (orb1146967)

Description	Human platelet-derived growth factor receptor, beta polypeptide (PDGFRB) is also called J03278, M216...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PDGF R beta
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 83.1 kDa. The protein migrates as 100-120 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_002600.1
Purity	95%
Source	Human PDGF R beta, Fc Tag (orb1146967) is expressed from human 293 cells (HEK293). It contains AA Leu 33 - Phe 530 (Accession # P09619-1).
MW	83.1 kDa
NCBI	NP_002600.1
Expiration Date	6 months from date of receipt.