


Product Datasheet

Human GFR alpha-like Protein (orb1146970)

Description	GFR alpha-like is also known as GDNF family receptor alpha-like, GFRAL, C6orf144, UNQ9356, PRO34128....
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	GFR alpha-like
Preservatives	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 64.2 kDa. The protein migrates as 70-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human GFR alpha-like, Fc Tag (orb1146970) is expressed from human 293 cells (HEK293). It contains AA Ser 19 - Glu 351 (Accession # Q6UXV0-1).
MW	64.2 kDa
NCBI	NP_997293.2
Expiration Date	6 months from date of receipt.



Immobilized
Human GDF-
15, His Tag
(Cat. ...)