

Product Datasheet

Mouse TNFR2 / CD120b / TNFRSF1B Protein (orb1146973)

Description

Mouse TNFR2, His Tag is expressed

contains AA Val 23 - Gly 258 (Ac...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TNFR2

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 27.2 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_035740.2

Purity

90%

Source

Mouse TNFR2, His Tag (orb1146973) is expressed from human 293 cells (HEK293). It contains AA Val 23 - Gly 258 (Accession # P25119-1).

MW

27.2 kDa

Uniprot ID

[P25119](#)

NCBI

[NP_035740.2](#)

Expiration Date

6 months from date of receipt.