

Product Datasheet

Mouse MIS RII Protein (orb1146988)

Description

Mouse MIS RII, His Tag is expressed from

Gln 16 - Met 147 (...)

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

MIS RII

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 15.9 kDa. The protein migrates as 22-24 kDa and 27-30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Mouse MIS RII, His Tag (orb1146988) is expressed from human 293 cells (HEK293). It contains AA Gln 16 - Met 147 (Accession # Q8K592-1).

MW

15.9 kDa

Uniprot ID

[Q8K592](#)

NCBI

[NP_653130.2](#)

Expiration Date

6 months from date of receipt.