

Product Datasheet

Human MIS RII Protein (orb1146989)

Description

Human MIS RII, His Tag is expressed from

Pro 18 - Ser 144 (...)

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

MIS RII

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 15.4 kDa. The protein migrates as 20-30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human MIS RII, His Tag (orb1146989) is expressed from human 293 cells (HEK293). It contains AA Pro 18 - Ser 144 (Accession # Q16671-1).

MW

15.4 kDa

Uniprot ID
Q16671
NCBI
NP_065434.1
Expiration Date

6 months from date of receipt.